Instructions for use

Philadelphia[®] Tracheotomy Collar

Life Without Limitations[®]

| | |
|-----------------------------------|----|
| | 3 |
| ⒺⒺ Instructions for use..... | 5 |
| ⒺⒺ Gebrauchsanweisung..... | 7 |
| ⒺⒺ Guide de fabrication..... | 9 |
| ⒺⒺ Instrucciones para el uso..... | 11 |
| ⒺⒺ Istruzioni per l'uso..... | 13 |
| ⒺⒺ Bruksanvisning..... | 15 |
| ⒺⒺ Brugsanvisning..... | 17 |
| ⒺⒺ Bruksanvisning..... | 19 |
| ⒺⒺ Käyttöohjeet..... | 21 |
| ⒺⒺ Gebruiksaanwijzing..... | 23 |
| ⒺⒺ Instruções de Utilização..... | 25 |

EN- This product has been designed and tested based on single patient usage. It is a single-use device and should never be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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
Measurement B - Neck Circumference

| Size | Circumference |
|-------------------|------------------------|
| Infant | 6"-8" (15.2-20.3 cm) |
| Pediatric | 8"-11" (20.3-27.9 cm) |
| Adult Small | 10"-13" (25.4-33.0 cm) |
| Adult Medium | 13"-16" (33.0-40.6 cm) |
| Adult Large | 16"-19" (40.6-48.3 cm) |
| Adult Extra Large | 19"+ (48.3 cm+) |

ENGLISH

MEASUREMENT & SIZING

Proper sizing is important for patient immobilization and comfort. Always maintain the patient's head in neutral alignment.

Height:

Measure the vertical distance from the top of the patient's shoulder to the tip of the chin (Figure A). The measurement corresponds to the HEIGHT of the collar. If the patient's measurement falls between two consecutive sizes, apply the smaller size first.

Circumference:

Measure around the patient's neck (Figure B). The measurement corresponds to the size CIRCUMFERENCE of the collar.

APPLICATION

Proper application of the collar is as important as proper sizing for immobilization and patient comfort. The two-piece design consists of a **FRONT** piece and a **BACK** piece which are packaged as a set.

1. After selecting the properly sized collar, apply the back piece of the collar to the back of the patient's neck. Center the collar. The back arrow should point upward.
2. Apply the **FRONT** piece of the collar with the chin secured in the recess. Center the collar to secure neutral alignment. The front piece **OVERLAPS** the back piece to ensure effective immobilization and comfort. The front arrow should point upward.
3. With hook and loop fasteners, tighten the collar with a bilateral adjustment. This will secure the patient's cervical region in neutral alignment.

FEATURES

- Non-toxic Plastazote® Foam
- Total Cervical Arch Support
- Full range of sizes for custom fit
- X-Ray Lucent
- Easily cleaned with soap and water
- Hook and loop fastener
- Use as prescribed by a physician

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