



THE JAS EZ TOE

ADVANCED TURNBUCKLE TECHNOLOGY

The JAS EZ Toe device offers full range, bi-directional ROM therapy for the great toe, from 90° dorsi-flexion to 45° plantar-flexion in a light weight, simple to use device. An innovative motion arc design, combined with infinitely adjustable ROM assures precise, pain-free, end range stretch for optimal therapeutic outcomes. Optional accessories are available to allow ROM therapy for 2nd through 5th toes.



JAS QUALITY AND PERFORMANCE IN A SINGLE-PATIENT USE DESIGN

Safe and Effective ROM Therapy Following:

- Bunionectomy
- Hallux Limitus
- Tendon and Ligament Repairs
- Toe Implants
- Arthroplasty
- Rheumatoid Arthritis

Infinite ROM Adjustment

Easy to apply, Simple to use

**JAS EZ Toe ROM
90° Dorsiflexion to
45° Plantarflexion**

Patented Motion Arm™ Eliminates Joint Compression

The JAS EZ Toe utilizes proven principles of stress relaxation to achieve permanent restoration of joint ROM with three 30-minute sessions per day.

Therapy is patient controlled with a simple turn of a knob, eliminating difficult screwdriver tension adjustments required with dynamic splints.

Features:	Benefits:
Patient controlled stretch via knob	Stretch adjusted without removing device or need for tools
Easy to apply and simple to use	Promotes excellent patient compliance
Works in both Dorsi and Plantar flexion	Full-range stress relaxation therapy in one device
Fulcrum positioned away from joint	Eliminates painful joint compression
Simulates manual stretch	Provides clinic-quality stretch in the home setting daily

AVAILABLE FOR COMMERCIAL PURCHASE



Joint Active Systems, Inc. / Telephone: 217.342.3412 or 800.879.0117 / Info@jointactivesystems.com

www.jointactivesystems.com

Licensed in the State of Illinois.