

A revolutionary way to reduce pain and swelling through the hot, cold, contrast and compression therapies you trust.

put a wrap on pain



VITALWEAR

VITALWRAP®  
SYSTEM

## A breakthrough approach to post-surgical pain: targeted relief, faster recoveries, better outcomes.

Put orthopedic patients on the fast track to recovery. The **SADMERC** - approved VitalWrap System combines hot, cold, contrast and compression therapies into one simple device. With our innovative family wraps, patients can treat pain and swelling with precision and comfort. VitalWrap's combined therapies work at every stage of rehabilitation, from the acute injury phase to post-surgery physical therapy. The result: better outcomes, faster recoveries and unprecedented flexibility.

### **SIMPLE:**

Portable, easy-to-use design for patient self-treatment.

### **PRECISE:**

Control temperature with one dial, 40° - 105° F.

### **FLEXIBLE:**

Hot, cold, contrast and compression therapies.

### **INNOVATIVE:**

Family of therapeutic wraps treats multiple body areas.



### **Wide-ranging applications — One versatile solution**

Arthritis • Back pain • Sprains, strains and inflammation • Post-surgical recovery  
Pain management • Functional mobility • Many other conditions



VITALWEAR

put a wrap on pain



### Systems

order no.

00002 VitalWrap System 120 V<sub>ac</sub>, 60 Hz  
(control unit, tubing set, standard VW)

00013 VitalWrap System 120 V<sub>ac</sub>, 60 Hz  
(control unit, tubing set, without standard VW)

### Wraps

00255	Medium VitalWrap (4.25" x 52")	for elbow and calf
00001	Standard VitalWrap (6" x 60")	general purpose wrap for calf, knee, hamstring, hip, back, chest and abdomen
00118	Large VitalWrap (8" x 72")	broad coverage wrap for leg and back
00388	Back VitalWrap	for use on back
00119	Medium Shoulder VitalWrap	for use on either shoulder
00395	Large Shoulder VitalWrap	for use on either shoulder
00394	Wrist VitalWrap	for hand and wrist
00396	Hand VitalWrap	for arthritic hands
00397	Medium Knee VitalWrap	for post-surgery and CPM applications
00392	Full Knee VitalWrap	for post-surgery and CPM applications
00398	Ankle VitalWrap	for foot and ankle
00122	Cervical VitalWrap	for neck applications
00648	Hip VitalWrap	for use on hip

### Accessories

00094	Carrying Case	for easy portability of the VW System
00034	VitalWrap Tubing Set	
00262	Hospital Grade Power Cord	

### SADMERC Durable Medical Equipment Codes

code

E0217	Water circulating heat pad with pump
E0218	Water circulating cold pad with pump
E0249	Pad for water circulation heat unit

## SYSTEM SPECIFICATIONS:

<b>PHYSICAL</b>	
size (approximate)	8" diameter x 13.5" height (20.3 cm x 34.3 cm)
weight (dry)	5 lbs (2.27 kg)
fluid tether	8' (2.44 m)
<b>CONTROL SYSTEM</b>	
type	manual user control
<b>THERMAL SYSTEM</b>	
range	40° - 105° F (4.4° - 40.6° C)
contrast therapy transition time	1 min typical from extremes in temperature range
operating period	2-6 hours typical continuous use
thermal cutout	116.6° F (47° C)
<b>CIRCULATING SYSTEM</b>	
reservoir capacity	1.1 gal (4.2 L)
reservoir fluid	ice water
flow rate thru VW	5.5 gph (21 Lph)
maximum pressure	18 psi
<b>ELECTRICAL SYSTEM:</b>	
voltage	120 V <sub>ac</sub> , 60 Hz
power	300 W
leakage current	under 300 µA
<b>CLASSIFICATION:</b>	
classification	class I equipment
type of equipment	type BF

00631 Rev. C

VitalWear, Inc.  
384 Oyster Point Blvd., Suite 16  
South San Francisco, CA 94080-1968  
Tel: 1-800-553-4081  
www.vitalwear.com