

THE JAS EZ KNEE EXTENSION

ADVANCED TURNBUCKLE TECHNOLOGY

The JAS EZ Knee Extension offers a no-compromise approach to Range of Motion (ROM) therapy – 48° flexion, 36° hyperextension – in a lightweight, easy-fit, single-patient use device. An innovative motion arc design, combined with infinitely adjustable ROM, assures precise end-range stretch throughout the entire treatment session. Maximal Total End Range Time (TERT) is achieved for optimal therapeutic benefit.



JAS QUALITY AND PERFORMANCE IN A SINGLE-PATIENT USE DESIGN

JAS EZ Knee Extension ROM
48° Flexion
36° Hyperextension

Optimal Comfort
and Compliance

Easy to Apply,
Simple to Use

Patented
Motion Arm™

The JAS EZ Knee Extension utilizes proven principles of stress relaxation to achieve permanent restoration of joint ROM in three 30-minute sessions per day.

The patient-controlled JAS protocol allows for *pain-free* therapy, virtually eliminating muscle guarding and assuring excellent compliance.

Features:	Benefits:
Anatomically correct arc of motion	Delivers optimal stretch loads through entire ROM
Lightweight, easy to apply and simple to use	Promotes excellent patient compliance
Fulcrum positioned away from the joint	Eliminates joint compression
Patient-controlled, infinitely adjustable ROM	Assures precise and pain-free end-range stretch; virtually eliminates muscle guarding
30-minute treatment sessions	Dramatically reduced treatment time ensures high patient compliance
Simulates manual stretch	Provides clinic-quality stretch in the home setting

AVAILABLE FOR COMMERCIAL PURCHASE



Joint Active Systems, Inc. / Telephone: 217.342.3412 or 800.879.0117 / Info@jointactivesystems.com

www.jointactivesystems.com

Licensed in the State of Illinois.

Mar 10



Safe and Effective ROM Therapy Following:

- Total Knee Replacement
- Ligament Repairs (ACL, PCL)
- Meniscectomy
- Post-Tibial/Femur Fractures
- Arthrofibrosis
- Burns

Knee Extension